

Endocervical Swab Specimen Collection Guide


Clinician-collection with the cobas[®] PCR Media Dual Swab Sample Kit

PROCEDURAL NOTES:

- Vaginal lubricants, speculum jellies, creams, and gels containing carbomer(s) may interfere with the test and should not be used during or prior to sample collection.
- If the collected specimen contains excess blood (specimen has a red or brown color), it should be discarded and not used for testing.
- Avoid contact of the cobas[®] PCR Media with the skin, eyes or mucous membranes. If contact does occur, immediately wash with large amounts of water.




NOTE: Do not pre-wet the swab in cobas[®] PCR Media before collection.



CLEAN:


- Using the woven swab, remove excess mucus from the cervical os and surrounding mucosa.
- Discard the swab after cleaning.

NOTE: Cleaning excess mucus from the cervical os is required to ensure an adequate sample is obtained for processing.




COLLECT:

- To collect the specimen, hold flocked swab with the scoreline above your hand and insert into the endocervical canal.
- Gently rotate the swab 5 times in one direction in the endocervical canal. Do not over-rotate.
- Carefully withdraw the swab, avoiding any contact with the vaginal mucosa.



ALIGN:

- Remove the cap from the cobas[®] PCR Media Tube.
- Lower the swab specimen into the tube until the visible scoreline on the swab shaft is aligned with the tube rim.
- The bud of the swab should not be submerged into the liquid prior to breaking the shaft.



BREAK:

- Carefully leverage the swab against the tube rim to break the swab shaft at the scoreline.

CLOSE:

- Tightly close the cobas[®] PCR Media Tube. The specimen is now ready for transport.
- Discard the top portion of the swab.

SPECIMEN TRANSPORT AND STORAGE:

- Transport and store the cobas[®] PCR Media Tube containing the collection swab at 2°C to 30°C.
- The specimen should only contain 1 flocked swab and may be rejected if the tube contains no swab or 2 swabs.