



Cardiac C-Reactive Protein (Latex) High Sensitive

cobas c 311/501/502 analyzers
04628918 190, Sys-ID 07 6866 9

For **cobas c 501/311 analyzers:**
CRPHS: ACN 217

For **cobas c 502 analyzer:**
CRPHS: ACN 8217

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
 Plasma / Plazma / Osocze / Πλάσμα

LOT	e-BC Ver.			mg/L			nmol/L			mg/dL		
	c311	c501	c502	Value	Range	1SD	Value	Range	1SD	Value	Range	1SD
				Wert	Bereich	1SD	Wert	Bereich	1SD	Wert	Bereich	1SD
				Valeur	Intervalle	1s	Valeur	Intervalle	1s	Valeur	Intervalle	1s
				Valor	Intervalo	1DE	Valor	Intervalo	1DE	Valor	Intervalo	1DE
				Valore	Intervallo	1DS	Valore	Intervallo	1DS	Valore	Intervallo	1DS
				Valor	Intervalo	1DP	Valor	Intervalo	1DP	Valor	Intervalo	1DP
				Värde	Område	1DP	Värde	Område	1DP	Värde	Område	1DP
				Værdi	Område	1SD	Værdi	Område	1SD	Værdi	Område	1SD
				Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD
				Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD
				Wartość	Zakres	1SD	Wartość	Zakres	1SD	Wartość	Zakres	1SD
				Érték	Tartomány	1SD	Érték	Tartomány	1SD	Érték	Tartomány	1SD
				Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

191102	2	2	2	9.31	7.90	-	10.72	0.47	88.6	75.4	-	101.8	4.4	0.931	0.790	-	1.072	0.047
191118	2	2	2	8.27	7.04	-	9.50	0.41	78.7	67.0	-	90.4	3.9	0.827	0.704	-	0.950	0.041

PreciControl ClinChem Multi 1, No. 05947626 160

191102	2	2	2	9.31	7.90	-	10.72	0.47	88.6	75.4	-	101.8	4.4	0.931	0.790	-	1.072	0.047
--------	---	---	---	------	------	---	-------	------	------	------	---	-------	-----	-------	-------	---	-------	-------