



Tina-quant C-Reactive Protein Gen.3

Roche/Hitachi analyzers

04956885 190

04956923 190

04995554 190

CRPL3: ACN 210 (MODULAR P)

CRPL3: ACN 200 (MODULAR D)

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
 Plasma / Plazma / Osocze / Πλάσμα

	mg/L	nmol/L	mg/dL
LOT	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s. Proteins, No. 11355279 216

144344	86.7	825	8.67
277667	86.7	825	8.67
199428	85.4	813	8.54
277668	85.4	813	8.54
281130	85.4	813	8.54
281134	85.4	813	8.54
281135	85.4	813	8.54
223235	80.9	770	8.09
277906	88.6	843	8.86
312505	88.0	838	8.80

C.f.a.s. Proteins, No. 11355279 160

144344	86.7	825	8.67
199428	85.4	813	8.54
223235	80.9	770	8.09
277906	88.6	843	8.86
312505	88.0	838	8.80

LOT	mg/L			nmol/L			mg/dL		
	Value Wert Valeur Valor Valore Valor Värde Værdi Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værdi Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værdi Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186380	8.63	6.89	-	10.37	0.58	82.2	65.7	-	98.7	5.5	0.863	0.689	-	1.037	0.058
132859	7.45	5.95	-	8.95	0.50	70.9	56.8	-	85.0	4.7	0.745	0.595	-	0.895	0.050
144514	7.37	5.90	-	8.84	0.49	70.2	56.1	-	84.3	4.7	0.737	0.590	-	0.884	0.049
154342	7.45	5.95	-	8.95	0.50	70.9	56.8	-	85.0	4.7	0.745	0.595	-	0.895	0.050
160400	7.54	6.04	-	9.04	0.50	71.8	57.4	-	86.2	4.8	0.754	0.604	-	0.904	0.050
160401	7.54	6.04	-	9.04	0.50	71.8	57.4	-	86.2	4.8	0.754	0.604	-	0.904	0.050
160407	8.41	6.73	-	10.09	0.56	80.1	64.2	-	96.0	5.3	0.841	0.673	-	1.009	0.056
160410	8.41	6.73	-	10.09	0.56	80.1	64.2	-	96.0	5.3	0.841	0.673	-	1.009	0.056
191102	8.20	6.55	-	9.85	0.55	78.1	62.5	-	93.7	5.2	0.820	0.655	-	0.985	0.055
191118	8.34	6.66	-	10.02	0.56	79.4	63.5	-	95.3	5.3	0.834	0.666	-	1.002	0.056
250241	8.76	7.02	-	10.50	0.58	83.4	66.6	-	100.2	5.6	0.876	0.702	-	1.050	0.058
250254	9.25	7.39	-	11.11	0.62	88.1	70.4	-	105.8	5.9	0.925	0.739	-	1.111	0.062

PreciControl ClinChem Multi 1, No. 05947626 160

144514	7.37	5.90	-	8.84	0.49	70.2	56.1	-	84.3	4.7	0.737	0.590	-	0.884	0.049
191102	8.20	6.55	-	9.85	0.55	78.1	62.5	-	93.7	5.2	0.820	0.655	-	0.985	0.055

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186394	37.5	30.0	-	45.0	2.5	357	285	-	429	24	3.75	3.00	-	4.50	0.25
186395	37.5	30.0	-	45.0	2.5	357	285	-	429	24	3.75	3.00	-	4.50	0.25
132865	36.2	29.0	-	43.4	2.4	345	276	-	414	23	3.62	2.90	-	4.34	0.24
144527	35.6	28.4	-	42.8	2.4	339	270	-	408	23	3.56	2.84	-	4.28	0.24
154345	36.3	29.1	-	43.5	2.4	346	277	-	415	23	3.63	2.91	-	4.35	0.24
160386	38.6	30.8	-	46.4	2.6	367	295	-	439	24	3.86	3.08	-	4.64	0.26
160393	39.0	31.2	-	46.8	2.6	371	296	-	446	25	3.90	3.12	-	4.68	0.26
191111	37.4	29.9	-	44.9	2.5	356	284	-	428	24	3.74	2.99	-	4.49	0.25
191112	37.4	29.9	-	44.9	2.5	356	284	-	428	24	3.74	2.99	-	4.49	0.25
191125	37.2	29.7	-	44.7	2.5	354	282	-	426	24	3.72	2.97	-	4.47	0.25
250228	53.9	43.1	-	64.7	3.6	513	411	-	615	34	5.39	4.31	-	6.47	0.36
250267	54.1	43.3	-	64.9	3.6	515	413	-	617	34	5.41	4.33	-	6.49	0.36

PreciControl ClinChem Multi 2, No. 05947774 160

186395	37.5	30.0	-	45.0	2.5	357	285	-	429	24	3.75	3.00	-	4.50	0.25
144527	35.6	28.4	-	42.8	2.4	339	270	-	408	23	3.56	2.84	-	4.28	0.24
154345	36.3	29.1	-	43.5	2.4	346	277	-	415	23	3.63	2.91	-	4.35	0.24

LOT	mg/L			nmol/L			mg/dL		
	Value	Range	1SD	Value	Range	1SD	Value	Range	1SD
	Wert	Bereich	1SD	Wert	Bereich	1SD	Wert	Bereich	1SD
	Valeur	Intervalle	1s	Valeur	Intervalle	1s	Valeur	Intervalle	1s
	Valor	Intervalo	1DE	Valor	Intervalo	1DE	Valor	Intervalo	1DE
	Valore	Intervallo	1DS	Valore	Intervallo	1DS	Valore	Intervallo	1DS
	Valor	Intervalo	1DP	Valor	Intervalo	1DP	Valor	Intervalo	1DP
	Värde	Område	1DP	Värde	Område	1DP	Värde	Område	1DP
	Værdi	Område	1SD	Værdi	Område	1SD	Værdi	Område	1SD
	Homota	Rozmezí	1SD	Homota	Rozmezí	1SD	Homota	Rozmezí	1SD
	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD
	Wartość	Zakres	1SD	Wartość	Zakres	1SD	Wartość	Zakres	1SD
	Érték	Tartomány	1SD	Érték	Tartomány	1SD	Érték	Tartomány	1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm Protein, No. 10557897 122

223237	9.69	7.74	-	11.64	0.65	92.2	73.9	-	110.5	6.1	0.969	0.774	-	1.164	0.065
294543	9.15	7.32	-	10.98	0.61	87.1	69.7	-	104.5	5.8	0.915	0.732	-	1.098	0.061

Precinorm Protein, No. 10557897 160

223237	9.69	7.74	-	11.64	0.65	92.2	73.9	-	110.5	6.1	0.969	0.774	-	1.164	0.065
294543	9.15	7.32	-	10.98	0.61	87.1	69.7	-	104.5	5.8	0.915	0.732	-	1.098	0.061

Precipath Protein, No. 11333127 122

154658	49.1	39.2	-	59.0	3.3	467	374	-	560	31	4.91	3.92	-	5.90	0.33
272990	50.3	40.1	-	60.5	3.4	479	383	-	575	32	5.03	4.01	-	6.05	0.34

Precipath Protein, No. 11333127 160

154658	49.1	39.2	-	59.0	3.3	467	374	-	560	31	4.91	3.92	-	5.90	0.33
272990	50.3	40.1	-	60.5	3.4	479	383	-	575	32	5.03	4.01	-	6.05	0.34