



C-Reactive Protein (Latex)

COBAS INTEGRA 400 plus analyzer

20764930 322, Sys.-ID 07 6493 0

CRPLX: Test-ID 0-293

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
 Plasma / Plazma / Osocze / Πλάσμα

	mg/L	nmol/L	mg/dL
LOT	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s. Proteins, No. 11355279 216

144344	104	990	10.4
277667	104	990	10.4
199428	99.6	948	9.96
277668	99.6	948	9.96
281130	99.6	948	9.96
281134	99.6	948	9.96
281135	99.6	948	9.96
223235	94.8	902	9.48
277906	105	1000	10.5
312505	105	1000	10.5

C.f.a.s. Proteins, No. 11355279 160

144344	104	990	10.4
199428	99.6	948	9.96
223235	94.8	902	9.48
277906	105	1000	10.5
312505	105	1000	10.5

LOT	mg/L			nmol/L			mg/dL		
	Value Wert Valeur Valor Valore Valor Värde Værði Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værði Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værði Hornota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD 1SD

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186380	9.16	7.33	-	10.99	0.61	87.2	69.8	-	104.6	5.8	0.916	0.733	-	1.099	0.061
132859	7.95	6.36	-	9.54	0.53	75.7	60.7	-	90.7	5.0	0.795	0.636	-	0.954	0.053
144514	7.74	6.18	-	9.30	0.52	73.7	59.0	-	88.4	4.9	0.774	0.618	-	0.930	0.052
154342	7.96	6.37	-	9.55	0.53	75.8	60.5	-	91.1	5.1	0.796	0.637	-	0.955	0.053
160400	8.00	6.41	-	9.59	0.53	76.2	60.9	-	91.5	5.1	0.800	0.641	-	0.959	0.053
160401	8.00	6.41	-	9.59	0.53	76.2	60.9	-	91.5	5.1	0.800	0.641	-	0.959	0.053
160407	8.87	7.10	-	10.64	0.59	84.4	67.6	-	101.2	5.6	0.887	0.710	-	1.064	0.059
160410	8.87	7.10	-	10.64	0.59	84.4	67.6	-	101.2	5.6	0.887	0.710	-	1.064	0.059
191102	8.88	7.11	-	10.65	0.59	84.5	67.7	-	101.3	5.6	0.888	0.711	-	1.065	0.059
191118	9.05	7.25	-	10.85	0.60	86.2	69.1	-	103.3	5.7	0.905	0.725	-	1.085	0.060
250241	9.54	7.62	-	11.46	0.64	90.8	72.5	-	109.1	6.1	0.954	0.762	-	1.146	0.064
250254	10.1	8.0	-	12.2	0.7	96.2	77.0	-	115.4	6.4	1.01	0.80	-	1.22	0.07

PreciControl ClinChem Multi 1, No. 05947626 160

144514	7.74	6.18	-	9.30	0.52	73.7	59.0	-	88.4	4.9	0.774	0.618	-	0.930	0.052
191102	8.88	7.11	-	10.65	0.59	84.5	67.7	-	101.3	5.6	0.888	0.711	-	1.065	0.059

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186394	39.0	31.2	-	46.8	2.6	371	296	-	446	25	3.90	3.12	-	4.68	0.26
186395	39.0	31.2	-	46.8	2.6	371	296	-	446	25	3.90	3.12	-	4.68	0.26
132865	39.0	31.2	-	46.8	2.6	371	296	-	446	25	3.90	3.12	-	4.68	0.26
144527	38.7	30.9	-	46.5	2.6	368	293	-	443	25	3.87	3.09	-	4.65	0.26
154345	38.8	31.0	-	46.6	2.6	369	294	-	444	25	3.88	3.10	-	4.66	0.26
160386	43.6	34.9	-	52.3	2.9	415	331	-	499	28	4.36	3.49	-	5.23	0.29
160393	43.1	34.4	-	51.8	2.9	410	329	-	491	27	4.31	3.44	-	5.18	0.29
191111	42.1	33.7	-	50.5	2.8	401	320	-	482	27	4.21	3.37	-	5.05	0.28
191112	42.1	33.7	-	50.5	2.8	401	320	-	482	27	4.21	3.37	-	5.05	0.28
191125	41.7	33.3	-	50.1	2.8	397	319	-	475	26	4.17	3.33	-	5.01	0.28
250228	60.1	48.1	-	72.1	4.0	572	458	-	686	38	6.01	4.81	-	7.21	0.40
250267	60.4	48.4	-	72.4	4.0	575	461	-	689	38	6.04	4.84	-	7.24	0.40

PreciControl ClinChem Multi 2, No. 05947774 160

186395	39.0	31.2	-	46.8	2.6	371	296	-	446	25	3.90	3.12	-	4.68	0.26
144527	38.7	30.9	-	46.5	2.6	368	293	-	443	25	3.87	3.09	-	4.65	0.26
154345	38.8	31.0	-	46.6	2.6	369	294	-	444	25	3.88	3.10	-	4.66	0.26

LOT	mg/L			nmol/L			mg/dL		
	Value	Range	1SD	Value	Range	1SD	Value	Range	1SD
	Wert	Bereich	1SD	Wert	Bereich	1SD	Wert	Bereich	1SD
	Valeur	Intervalle	1s	Valeur	Intervalle	1s	Valeur	Intervalle	1s
	Valor	Intervalo	1DE	Valor	Intervalo	1DE	Valor	Intervalo	1DE
	Valore	Intervallo	1DS	Valore	Intervallo	1DS	Valore	Intervallo	1DS
	Valor	Intervalo	1DP	Valor	Intervalo	1DP	Valor	Intervalo	1DP
	Värde	Område	1DP	Värde	Område	1DP	Värde	Område	1DP
	Værdi	Område	1SD	Værdi	Område	1SD	Værdi	Område	1SD
	Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD
	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD
	Wartość	Zakres	1SD	Wartość	Zakres	1SD	Wartość	Zakres	1SD
	Érték	Tartomány	1SD	Érték	Tartomány	1SD	Érték	Tartomány	1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm Protein, No. 10557897 122

223237	12.1	9.7	-	14.5	0.8	115	91	-	139	8	1.21	0.97	-	1.45	0.08
294543	11.4	9.0	-	13.8	0.8	109	88	-	130	7	1.14	0.90	-	1.38	0.08

Precinorm Protein, No. 10557897 160

223237	12.1	9.7	-	14.5	0.8	115	91	-	139	8	1.21	0.97	-	1.45	0.08
294543	11.4	9.0	-	13.8	0.8	109	88	-	130	7	1.14	0.90	-	1.38	0.08

Precipath Protein, No. 11333127 122

154658	60.8	48.5	-	73.1	4.1	579	462	-	696	39	6.08	4.85	-	7.31	0.41
272990	60.4	48.4	-	72.4	4.0	575	461	-	689	38	6.04	4.84	-	7.24	0.40

Precipath Protein, No. 11333127 160

154658	60.8	48.5	-	73.1	4.1	579	462	-	696	39	6.08	4.85	-	7.31	0.41
272990	60.4	48.4	-	72.4	4.0	575	461	-	689	38	6.04	4.84	-	7.24	0.40