



C-reactive protein (latex)

cobas c 111 analyzer

05401615 190

CRPLX: ACN 019

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
 Plasma / Plazma / Osocze / Πλάσμα

	mg/L	nmol/L
LOT	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s. Proteins, No. 11355279 216

144344	104	990
277667	104	990
199428	99.6	948
277668	99.6	948
281130	99.6	948
281134	99.6	948
281135	99.6	948
223235	94.8	902
277906	105	1000
312505	105	1000

C.f.a.s. Proteins, No. 11355279 160

144344	104	990
199428	99.6	948
223235	94.8	902
277906	105	1000
312505	105	1000

LOT	mg/L			nmol/L		
	Value	Range	1SD	Value	Range	1SD
Wert	Bereich	1SD	Wert	Bereich	1SD	
Valeur	Intervalle	1s	Valeur	Intervalle	1s	
Valor	Intervalo	1DE	Valor	Intervalo	1DE	
Valore	Intervallo	1DS	Valore	Intervallo	1DS	
Valor	Intervalo	1DP	Valor	Intervalo	1DP	
Värde	Område	1DP	Värde	Område	1DP	
Værdi	Område	1SD	Værdi	Område	1SD	
Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD	
Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	
Wartość	Zakres	1SD	Wartość	Zakres	1SD	
Érték	Tartomány	1SD	Érték	Tartomány	1SD	
Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186380	9.54	7.62	-	11.46	0.64	90.8	72.5	-	109.1	6.1
132859	8.28	6.63	-	9.93	0.55	78.8	62.9	-	94.7	5.3
144514	8.06	6.44	-	9.68	0.54	76.7	61.4	-	92.0	5.1
154342	8.29	6.64	-	9.94	0.55	78.9	63.0	-	94.8	5.3
160400	8.33	6.65	-	10.01	0.56	79.3	63.4	-	95.2	5.3
160401	8.33	6.65	-	10.01	0.56	79.3	63.4	-	95.2	5.3
160407	9.23	7.37	-	11.09	0.62	87.9	70.2	-	105.6	5.9
160410	9.23	7.37	-	11.09	0.62	87.9	70.2	-	105.6	5.9
191102	9.24	7.38	-	11.10	0.62	88.0	70.3	-	105.7	5.9
191118	9.42	7.53	-	11.31	0.63	89.7	71.7	-	107.7	6.0
250241	9.93	7.95	-	11.91	0.66	94.5	75.6	-	113.4	6.3
250254	10.5	8.4	-	12.6	0.7	100	79	-	121	7

PreciControl ClinChem Multi 1, No. 05947626 160

144514	8.06	6.44	-	9.68	0.54	76.7	61.4	-	92.0	5.1
191102	9.24	7.38	-	11.10	0.62	88.0	70.3	-	105.7	5.9

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186394	40.0	31.9	-	48.1	2.7	381	306	-	456	25
186395	40.0	31.9	-	48.1	2.7	381	306	-	456	25
132865	40.0	31.9	-	48.1	2.7	381	306	-	456	25
144527	39.7	31.9	-	47.5	2.6	378	303	-	453	25
154345	39.8	31.7	-	47.9	2.7	379	304	-	454	25
160386	44.7	35.7	-	53.7	3.0	426	342	-	510	28
160393	44.2	35.5	-	52.9	2.9	421	337	-	505	28
191111	43.2	34.5	-	51.9	2.9	411	330	-	492	27
191112	43.2	34.5	-	51.9	2.9	411	330	-	492	27
191125	42.7	34.3	-	51.1	2.8	407	326	-	488	27
250228	61.6	49.3	-	73.9	4.1	586	469	-	703	39
250267	61.9	49.6	-	74.2	4.1	589	472	-	706	39

PreciControl ClinChem Multi 2, No. 05947774 160

186395	40.0	31.9	-	48.1	2.7	381	306	-	456	25
144527	39.7	31.9	-	47.5	2.6	378	303	-	453	25
154345	39.8	31.7	-	47.9	2.7	379	304	-	454	25

CRPLX

cobas c 111 analyzer

cobas®

Note 2018-06-05A01615 - 3 - 01

LOT	mg/L			nmol/L		
	Value	Range	1SD	Value	Range	1SD
Wert	Bereich	1SD	Wert	Bereich	1SD	
Valeur	Intervalle	1s	Valeur	Intervalle	1s	
Valor	Intervalo	1DE	Valor	Intervalo	1DE	
Valore	Intervallo	1DS	Valore	Intervallo	1DS	
Valor	Intervalo	1DP	Valor	Intervalo	1DP	
Värde	Område	1DP	Värde	Område	1DP	
Værdi	Område	1SD	Værdi	Område	1SD	
Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD	
Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	
Wartość	Zakres	1SD	Wartość	Zakres	1SD	
Érték	Tartomány	1SD	Érték	Tartomány	1SD	
Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	

Precinorm Protein, No. 10557897 122

223237	12.7	10.3	-	15.1	0.8	121	97	-	145	8
294543	11.9	9.5	-	14.3	0.8	113	89	-	137	8

Precinorm Protein, No. 10557897 160

223237	12.7	10.3	-	15.1	0.8	121	97	-	145	8
294543	11.9	9.5	-	14.3	0.8	113	89	-	137	8

Precipath Protein, No. 11333127 122

154658	63.2	50.6	-	75.8	4.2	602	482	-	722	40
272990	62.8	50.2	-	75.4	4.2	598	478	-	718	40

Precipath Protein, No. 11333127 160

154658	63.2	50.6	-	75.8	4.2	602	482	-	722	40
272990	62.8	50.2	-	75.4	4.2	598	478	-	718	40

2018-06

