



Tina-quant IgG Gen.2

cobas c 311 analyzer
03507432 190, Sys-ID 07 6787 5

cobas c 701/702 analyzers
05220718 190, Sys-ID 03 6787 5

For **cobas c 311 analyzer**:
IGG-2: ACN 674

For **cobas c 701/702 analyzers**:
IGG-2: ACN 8674

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
Plasma / Plazma / Osocze / Πλάσμα

LOT	e-BC Ver.			g/L	μmol/L	mg/dL
	c311	c701	c702	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsværdi/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsværdi/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsværdi/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

C.f.a.s. Proteins, No. 11355279 216

144344	3	3	3	19.1	127	1910
277667	2	2	2	19.1	127	1910
199428	3	3	3	18.8	125	1880
277668	2	2	2	18.8	125	1880
281130	2	2	2	18.8	125	1880
281134	2	2	2	18.8	125	1880
281135	2	2	2	18.8	125	1880
223235	2	2	2	19.7	131	1970
277906	2	2	2	19.5	130	1950
312505	1	1	1	19.6	131	1960

C.f.a.s. Proteins, No. 11355279 160

144344	3	3	3	19.1	127	1910
199428	3	3	3	18.8	125	1880
223235	2	2	2	19.7	131	1970
277906	2	2	2	19.5	130	1950
312505	1	1	1	19.6	131	1960

LOT	e-BC Ver.			g/L			μmol/L			mg/dL		
	c311	c701	c702	Value Wert Valeur Valor Valore Valor Värde Værdi Homota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værdi Homota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD	Value Wert Valeur Valor Valore Valor Värde Værdi Homota Hodnota Wartość Érték Τιμή	Range Bereich Intervalle Intervalo Intervallo Intervalo Område Område Rozmezí Rozsah Zakres Tartomány Εύρος	1SD 1SD 1s 1DE 1DS 1DP 1DP 1SD 1SD 1SD 1SD 1SD

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186380	5	5	5	7.40	6.08 - 8.72	0.44	49.4	40.4 - 58.4	3.0	740	608 - 872	44
132859	3	3	3	7.28	5.96 - 8.60	0.44	48.6	39.9 - 57.3	2.9	728	596 - 860	44
144514	4	4	4	7.23	5.94 - 8.52	0.43	48.2	39.5 - 56.9	2.9	723	594 - 852	43
154342	4	4	4	7.27	5.95 - 8.59	0.44	48.5	39.8 - 57.2	2.9	727	595 - 859	44
160400	3	3	3	7.12	5.83 - 8.41	0.43	47.5	38.8 - 56.2	2.9	712	583 - 841	43
160401	3	3	3	7.12	5.83 - 8.41	0.43	47.5	38.8 - 56.2	2.9	712	583 - 841	43
160407	4	4	4	7.31	5.99 - 8.63	0.44	48.8	40.1 - 57.5	2.9	731	599 - 863	44
160410	4	4	4	7.31	5.99 - 8.63	0.44	48.8	40.1 - 57.5	2.9	731	599 - 863	44
191102	3	2	2	6.97	5.71 - 8.23	0.42	46.5	38.1 - 54.9	2.8	697	571 - 823	42
191118	3	2	2	7.02	5.76 - 8.28	0.42	46.8	38.4 - 55.2	2.8	702	576 - 828	42
250241	3	3	3	7.23	5.94 - 8.52	0.43	48.2	39.5 - 56.9	2.9	723	594 - 852	43
250254	1	1	1	7.49	6.14 - 8.84	0.45	50.0	41.0 - 59.0	3.0	749	614 - 884	45

PreciControl ClinChem Multi 1, No. 05947626 160

144514	4	4	4	7.23	5.94 - 8.52	0.43	48.2	39.5 - 56.9	2.9	723	594 - 852	43
191102	3	2	2	6.97	5.71 - 8.23	0.42	46.5	38.1 - 54.9	2.8	697	571 - 823	42

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186394	6	6	6	12.0	9.9 - 14.1	0.7	80.0	65.6 - 94.4	4.8	1200	990 - 1410	70
186395	6	6	6	12.0	9.9 - 14.1	0.7	80.0	65.6 - 94.4	4.8	1200	990 - 1410	70
132865	5	5	5	11.7	9.6 - 13.8	0.7	78.0	63.9 - 92.1	4.7	1170	960 - 1380	70
144527	5	5	5	11.2	9.1 - 13.3	0.7	74.7	61.2 - 88.2	4.5	1120	910 - 1330	70
154345	5	4	4	11.7	9.6 - 13.8	0.7	78.0	63.9 - 92.1	4.7	1170	960 - 1380	70
160386	4	4	4	10.7	8.9 - 12.5	0.6	71.4	58.5 - 84.3	4.3	1070	890 - 1250	60
160393	5	5	5	10.9	8.8 - 13.0	0.7	72.7	59.5 - 85.9	4.4	1090	880 - 1300	70
191111	3	3	3	11.0	8.9 - 13.1	0.7	73.4	60.2 - 86.6	4.4	1100	890 - 1310	70
191112	3	3	3	11.0	8.9 - 13.1	0.7	73.4	60.2 - 86.6	4.4	1100	890 - 1310	70
191125	3	3	3	10.6	8.8 - 12.4	0.6	70.7	58.1 - 83.3	4.2	1060	880 - 1240	60
250228	3	3	3	12.0	9.9 - 14.1	0.7	80.0	65.6 - 94.4	4.8	1200	990 - 1410	70
250267	1	1	1	11.5	9.4 - 13.6	0.7	76.7	62.9 - 90.5	4.6	1150	940 - 1360	70

PreciControl ClinChem Multi 2, No. 05947774 160

186395	6	6	6	12.0	9.9 - 14.1	0.7	80.0	65.6 - 94.4	4.8	1200	990 - 1410	70
144527	5	5	5	11.2	9.1 - 13.3	0.7	74.7	61.2 - 88.2	4.5	1120	910 - 1330	70
154345	3	3	3	11.7	9.6 - 13.8	0.7	78.0	63.9 - 92.1	4.7	1170	960 - 1380	70

LOT	e-BC Ver.			g/L			μmol/L			mg/dL		
	c311	c701	c702	Value	Range	1SD	Value	Range	1SD	Value	Range	1SD
				Wert	Bereich	1SD	Wert	Bereich	1SD	Wert	Bereich	1SD
				Valeur	Intervalle	1s	Valeur	Intervalle	1s	Valeur	Intervalle	1s
				Valor	Intervalo	1DE	Valor	Intervalo	1DE	Valor	Intervalo	1DE
				Valore	Intervallo	1DS	Valore	Intervallo	1DS	Valore	Intervallo	1DS
				Valor	Intervalo	1DP	Valor	Intervalo	1DP	Valor	Intervalo	1DP
				Vårde	Område	1DP	Vårde	Område	1DP	Vårde	Område	1DP
				Værdi	Område	1SD	Værdi	Område	1SD	Værdi	Område	1SD
				Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD	Hornota	Rozmezí	1SD
				Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD	Hodnota	Rozsah	1SD
				Wartość	Zakres	1SD	Wartość	Zakres	1SD	Wartość	Zakres	1SD
				Érték	Tartomány	1SD	Érték	Tartomány	1SD	Érték	Tartomány	1SD
				Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm Protein, No. 10557897 122

223237	2	2	2	8.59	7.03 - 10.15	0.52	57.3	47.1 - 67.5	3.4	859	703 - 1015	52
294543	1	1	1	8.98	7.36 - 10.60	0.54	59.9	49.1 - 70.7	3.6	898	736 - 1060	54

Precinorm Protein, No. 10557897 160

223237	2	2	2	8.59	7.03 - 10.15	0.52	57.3	47.1 - 67.5	3.4	859	703 - 1015	52
294543	1	1	1	8.98	7.36 - 10.60	0.54	59.9	49.1 - 70.7	3.6	898	736 - 1060	54

Precipath Protein, No. 11333127 122

154658	2	2	2	16.6	13.6 - 19.6	1.0	111	90 - 132	7	1660	1360 - 1960	100
272990	1	1	1	17.4	14.4 - 20.4	1.0	116	95 - 137	7	1740	1440 - 2040	100

Precipath Protein, No. 11333127 160

154658	2	2	2	16.6	13.6 - 19.6	1.0	111	90 - 132	7	1660	1360 - 1960	100
272990	1	1	1	17.4	14.4 - 20.4	1.0	116	95 - 137	7	1740	1440 - 2040	100