

Supplemental Value Sheet



COBAS INTEGRA 400 plus analyzer

Please note that following data have been modified or added.

CERU Test-ID 0-666 Page 2-3

CERU



Ceruloplasmin

COBAS INTEGRA 400 plus analyzer


20764663 322, System-ID 07 6466 3

CERU3: Test-ID 0-666

LOT		g/L	$\mu\text{mol/L}$	mg/dL
		Calibration Value		Calibration Value

C.f.a.s. PAC, No. 03555941 190

305634	2020-07	0.333	2.48	33.3
383319	2021-05	0.318	2.37	31.8

LOT		g/L			µmol/L			mg/dL		
		Value	Range	1SD	Value	Range	1SD	Value	Range	1SD

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

191118	2020-03	0.134	0.101 - 0.167	0.011	1.00	0.76 - 1.24	0.08	13.4	10.1 - 16.7	1.1
250241	2020-11	0.136	0.103 - 0.169	0.011	1.01	0.77 - 1.25	0.08	13.6	10.3 - 16.9	1.1
250254	2020-10	0.145	0.109 - 0.181	0.012	1.08	0.81 - 1.35	0.09	14.5	10.9 - 18.1	1.2
324196	2021-04	0.152	0.116 - 0.188	0.012	1.13	0.86 - 1.40	0.09	15.2	11.6 - 18.8	1.2
324209	2021-08	0.156	0.120 - 0.192	0.012	1.16	0.89 - 1.43	0.09	15.6	12.0 - 19.2	1.2
324213	2021-07	0.160	0.121 - 0.199	0.013	1.19	0.89 - 1.49	0.10	16.0	12.1 - 19.9	1.3
324222	2021-10	0.138	0.105 - 0.171	0.011	1.03	0.79 - 1.27	0.08	13.8	10.5 - 17.1	1.1
410037	2022-07	0.151	0.115 - 0.187	0.012	1.13	0.86 - 1.40	0.09	15.1	11.5 - 18.7	1.2

PreciControl ClinChem Multi 1, No. 05947626 160

250241	2020-11	0.136	0.103 - 0.169	0.011	1.01	0.77 - 1.25	0.08	13.6	10.3 - 16.9	1.1
250254	2020-10	0.145	0.109 - 0.181	0.012	1.08	0.81 - 1.35	0.09	14.5	10.9 - 18.1	1.2
324209	2021-08	0.156	0.120 - 0.192	0.012	1.16	0.89 - 1.43	0.09	15.6	12.0 - 19.2	1.2
410037	2022-07	0.151	0.115 - 0.187	0.012	1.13	0.86 - 1.40	0.09	15.1	11.5 - 18.7	1.2

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

191125	2020-03	0.234	0.177 - 0.291	0.019	1.75	1.33 - 2.17	0.14	23.4	17.7 - 29.1	1.9
250228	2020-06	0.266	0.203 - 0.329	0.021	1.98	1.50 - 2.46	0.16	26.6	20.3 - 32.9	2.1
250267	2020-11	0.247	0.187 - 0.307	0.020	1.84	1.39 - 2.29	0.15	24.7	18.7 - 30.7	2.0
250280	2020-12	0.206	0.158 - 0.254	0.016	1.54	1.18 - 1.90	0.12	20.6	15.8 - 25.4	1.6
324345	2021-04	0.220	0.166 - 0.274	0.018	1.64	1.25 - 2.03	0.13	22.0	16.6 - 27.4	1.8
324348	2021-04	0.220	0.166 - 0.274	0.018	1.64	1.25 - 2.03	0.13	22.0	16.6 - 27.4	1.8
410044	2022-07	0.257	0.194 - 0.320	0.021	1.92	1.47 - 2.37	0.15	25.7	19.4 - 32.0	2.1
410045	2022-07	0.258	0.195 - 0.321	0.021	1.92	1.47 - 2.37	0.15	25.8	19.5 - 32.1	2.1

PreciControl ClinChem Multi 2, No. 05947774 160

250267	2020-11	0.247	0.187 - 0.307	0.020	1.84	1.39 - 2.29	0.15	24.7	18.7 - 30.7	2.0
250280	2020-12	0.206	0.158 - 0.254	0.016	1.54	1.18 - 1.90	0.12	20.6	15.8 - 25.4	1.6
324345	2021-04	0.220	0.166 - 0.274	0.018	1.64	1.25 - 2.03	0.13	22.0	16.6 - 27.4	1.8
324348	2021-04	0.220	0.166 - 0.274	0.018	1.64	1.25 - 2.03	0.13	22.0	16.6 - 27.4	1.8

Supplemental Value Sheet

Glossary

EN	Glossary	Calibration Value	Value	Range	Supplemental Value Sheet
DE	Glossar	Kalibrationswert	Wert	Bereich	Ergänzendes Wertebblatt
FR	Glossaire	Valeur de calibration	Valeur	Intervalle	Fiche de valeurs complémentaire
ES	Glosario	Valor de calibración	Valor	Intervalo	Hoja adicional de valores
IT	Glossario	Valore di calibrazione	Valore	Intervallo	Foglietto aggiuntivo relativo ai valori teorici
PT	Glossário	Valor de calibração	Valor	Intervalo	Folha de Valores Suplementar
DA	Ordliste	Kalibreringsværdi	Værdi	Område	Supplerende værdiskema
SV	Ordlista	Kalibreringsvärde	Värde	Intervall	Kompletterande värdeark
NO	Ordliste	Kalibreringsverdi	Verdi	Område	supplerende verdiark
CS	Glosář	Kalibrační hodnota	Hodnota	Rozmezí	Doplňkový list hodnot
SK	Slovník	Kalibračná hodnota	Hodnota	Rozsah	Doplňková tabuľka hodnôt
PL	Słownik	Wartość kalibracji	Wartość	Zakres	Uzupełniający arkusz wartości
HU	Fogalomtár	Kalibrációs érték	Érték	Tartomány	pótadatlap
EL	Γλωσσάριο	Τιμή βαθμονόμησης	Τιμή	Εύρος	Συμπληρωματικό φύλλο τιμών
TR	Sözlükçe	Kalibrasyon değeri	Değer	Aralık	Tamamlayıcı Değer Sayfası
BG	Речник	Калибрационна стойност	Стойност	Обхват	Допълнителна листовка със стойности
ET	Sõnastik	Kalibreerimisväärtus	Väärtus	Vahemik	Lisaväärtuste leht
LT	Terminų žodynas	Kalibravimo reikšmė	Reikšmė	Intervalas	Papildoma reikšmių lentelė
LV	Vārdnīca	Kalibrācijas vērtība	Vērtība	Diapazons	Papildu vērtību lapa
RU	Глоссарий	Калибровочное значение	Значение	Диапазон	Дополнительный перечень значений
VI	Bảng chú giải thuật ngữ	Giá trị chuẩn định	Giá trị	Khoảng	Tờ giá trị bổ sung
ZH	词汇表	定标值	数值	范围	补充靶值单

EN Please note that following data have been modified or added.

DE Bitte beachten Sie, dass folgende Daten geändert oder hinzugefügt wurden.

FR Veuillez noter que les données suivantes ont été modifiées ou ajoutées.

ES Observe que los siguientes datos han sido modificados o añadidos.

IT Nota: i seguenti dati sono stati modificati o aggiunti.

PT Note que os seguintes dados foram modificados ou adicionados.

DA Bemærk, at nedenstående data er ændret eller tilføjet.

SV Observera att följande data har ändrats eller lagts till.

NO Vær oppmerksom på at følgende data er endret eller lagt til.

CS Vezměte, prosím, na vědomí, že následující údaje byly změněny nebo přidány.

SK Upozorňujeme Vás, že nasledujúce údaje boli upravené alebo pridané.

PL Należy wziąć pod uwagę, że podane poniżej dane zostały zmodyfikowane lub zostały dodane.

HU Kérjük, vegyék figyelembe, hogy az alábbi adatokat módosították ill. hozzáadták.

EL Λάβετε υπόψη ότι τα παρακάτω δεδομένα έχουν τροποποιηθεί ή προστεθεί.

TR Aşağıdaki verilerde değişiklik veya ekleme yapıldığını lütfen dikkate alın.

BG Моля, обърнете внимание, че следните данни са променени или добавени.

ET Juhime tähelepanu, et järgmised andmed on muudetud või lisatud.

LT Atkreipkite dėmesį, kad šie duomenys buvo pakeisti arba papildyti.

LV Lūdzu, ņemiet vērā, ka ir mainīti vai pievienoti turpmāk norādītie dati.

RU Обратите внимание, что следующие данные были изменены или добавлены.

VI Vui lòng lưu ý rằng dữ liệu sau đã được sửa đổi hoặc bổ sung.

ZH 请注意：已修改或添加下列数据。

2020-02

