



Lipase colorimetric assay

cobas c 311/501/502 analyzers
03029590 322, Sys-ID 07 5900 7

cobas c 701/702 analyzers
07041918 190, Sys-ID 05 5900 8

For **cobas c 501/311 analyzers**:
LIPC: ACN 789 LIPC: ACN 786 (STAT)

For **cobas c 502/701/702 analyzers**:
LIPC: ACN 8789 LIPC: ACN 8786 (STAT)

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
Plasma / Plazma / Osocze / Πλάσμα

LOT	e-BC Ver.					U/L	µkat/L
	c311	c501	c502	c701	c702	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärði Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärði Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

Cfas, No. 10759350 190

186422	2	2	2	2	2	98.9	1.65
186423	2	2	2	2	2	99.5	1.66
186424	3	3	3	3	3	98.5	1.64
186425	3	3	3	3	3	98.6	1.65
186426	3	3	3	3	3	98.5	1.64
186427	3	3	3	3	3	99.5	1.66
186487	3	3	3	3	3	102	1.70
186488	3	3	3	3	3	102	1.70
186489	2	2	2	2	2	102	1.70
186490	3	3	3	3	3	102	1.70
144455	2	2	2	2	2	104	1.74
144456	2	2	2	2	2	105	1.75
144457	2	2	2	2	2	104	1.74
198885	2	2	2	2	2	107	1.79
198886	2	2	2	2	2	108	1.80
198888	2	2	2	2	2	107	1.79
198889	2	2	2	2	2	107	1.79
198890	2	2	2	2	2	107	1.79

LOT	e-BC Ver.					U/L			µkat/L		
	c311	c501	c502	c701	c702	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

Precinorm U, No. 10171735 122, 10171743 122

176367	4	4	4	4	4	45.4	37.3 - 53.5	2.7	0.758	0.623 - 0.893	0.045
186436	4	4	4	4	4	44.6	36.5 - 52.7	2.7	0.745	0.610 - 0.880	0.045
186438	2	2	2	2	2	44.6	36.5 - 52.7	2.7	0.745	0.610 - 0.880	0.045
186439	2	2	2	2	2	44.6	36.5 - 52.7	2.7	0.745	0.610 - 0.880	0.045

Precinorm U plus, No. 12149435 122

186350	2	2	2	2	2	43.8	36.0 - 51.6	2.6	0.731	0.599 - 0.863	0.044
203067	1	1	1	1	1	45.5	37.4 - 53.6	2.7	0.760	0.622 - 0.898	0.046
203068	1	1	1	1	1	45.3	37.2 - 53.4	2.7	0.757	0.622 - 0.892	0.045

Precipath U, No. 10171760 122, 10171778 122

176371	4	4	4	4	4	104	86 - 122	6	1.74	1.44 - 2.04	0.10
186443	4	4	4	4	4	108	90 - 126	6	1.80	1.47 - 2.13	0.11

Precipath U plus, No. 12149443 122

186345	2	2	2	2	2	94.9	77.8 - 112.0	5.7	1.58	1.31 - 1.85	0.09
203060	2	2	2	2	2	99.2	81.2 - 117.2	6.0	1.66	1.36 - 1.96	0.10

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186373	5	5	5	5	5	43.7	35.9 - 51.5	2.6	0.730	0.598 - 0.862	0.044
186380	4	4	4	4	4	49.0	40.3 - 57.7	2.9	0.818	0.671 - 0.965	0.049
132859	2	2	2	2	2	44.6	36.5 - 52.7	2.7	0.745	0.610 - 0.880	0.045
154342	3	3	3	3	3	44.8	36.7 - 52.9	2.7	0.748	0.613 - 0.883	0.045
144514	3	3	3	3	3	45.8	37.7 - 53.9	2.7	0.765	0.627 - 0.903	0.046
160400	2	2	2	2	2	45.6	37.5 - 53.7	2.7	0.762	0.624 - 0.900	0.046
160401	2	2	2	2	2	45.6	37.5 - 53.7	2.7	0.762	0.624 - 0.900	0.046
160407	3	3	3	3	3	47.4	39.0 - 55.8	2.8	0.792	0.648 - 0.936	0.048
160410	3	3	3	3	3	47.4	39.0 - 55.8	2.8	0.792	0.648 - 0.936	0.048
191102	1	1	1	1	1	45.3	37.2 - 53.4	2.7	0.757	0.622 - 0.892	0.045
191118	1	1	1	1	1	47.3	38.9 - 55.7	2.8	0.790	0.649 - 0.931	0.047

LOT	e-BC Ver.					U/L			µkat/L		
	c311	c501	c502	c701	c702	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186387	5	5	5	5	5	98.3	80.6 - 116.0	5.9	1.64	1.34 - 1.94	0.10
186394	4	4	4	4	4	100	82 - 118	6	1.67	1.37 - 1.97	0.10
186395	4	4	4	4	4	100	82 - 118	6	1.67	1.37 - 1.97	0.10
154345	3	2	2	2	2	99.1	81.4 - 116.8	5.9	1.65	1.35 - 1.95	0.10
144527	3	3	3	3	3	100	82 - 118	6	1.67	1.37 - 1.97	0.10
132865	3	3	3	3	3	101	83 - 119	6	1.69	1.39 - 1.99	0.10
160386	2	2	2	2	2	99.1	81.4 - 116.8	5.9	1.65	1.35 - 1.95	0.10
160393	3	3	3	3	3	102	84 - 120	6	1.70	1.40 - 2.00	0.10
191111	1	1	1	1	1	98.7	81.0 - 116.4	5.9	1.65	1.35 - 1.95	0.10
191112	1	1	1	1	1	98.7	81.0 - 116.4	5.9	1.65	1.35 - 1.95	0.10
191125	1	1	1	1	1	101	83 - 119	6	1.69	1.39 - 1.99	0.10