

COBAS INTEGRA 400 plus/800

Note 2017-06-03333752 - 2-3

LOT	U/L			µkat/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm U, No. 10171735 122, 10171743 122

176367	84.8	69.5 - 100.1	5.1	1.42	1.15 - 1.69	0.09
186436	89.9	73.7 - 106.1	5.4	1.50	1.23 - 1.77	0.09
186438	95.0	77.9 - 112.1	5.7	1.59	1.29 - 1.89	0.10
186439	94.1	77.3 - 110.9	5.6	1.57	1.30 - 1.84	0.09

Precinorm U plus, No. 12149435 122

186350	85.5	70.2 - 100.8	5.1	1.43	1.16 - 1.70	0.09
203067	80.1	65.7 - 94.5	4.8	1.34	1.10 - 1.58	0.08
203068	79.4	65.0 - 93.8	4.8	1.33	1.09 - 1.57	0.08

Precipath U, No. 10171760 122, 10171778 122

176371	228	186 - 270	14	3.81	3.12 - 4.50	0.23
186443	237	195 - 279	14	3.96	3.24 - 4.68	0.24

Precipath U plus, No. 12149443 122

186345	229	187 - 271	14	3.82	3.13 - 4.51	0.23
203060	221	182 - 260	13	3.69	3.03 - 4.35	0.22

PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)

186373	95.9	78.5 - 113.3	5.8	1.60	1.30 - 1.90	0.10
186380	94.1	77.3 - 110.9	5.6	1.57	1.30 - 1.84	0.09
132859	96.7	79.3 - 114.1	5.8	1.61	1.31 - 1.91	0.10
154342	93.2	76.4 - 110.0	5.6	1.56	1.29 - 1.83	0.09
144514	91.7	75.2 - 108.2	5.5	1.53	1.26 - 1.80	0.09
160400	98.9	81.2 - 116.6	5.9	1.65	1.35 - 1.95	0.10
160401	100	82 - 118	6	1.67	1.37 - 1.97	0.10
160407	106	88 - 124	6	1.77	1.44 - 2.10	0.11
160410	105	87 - 123	6	1.75	1.42 - 2.08	0.11
191102	93.6	76.8 - 110.4	5.6	1.56	1.29 - 1.83	0.09
191118	107	89 - 125	6	1.79	1.46 - 2.12	0.11

PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)

186387	229	187 - 271	14	3.82	3.13 - 4.51	0.23
186394	204	168 - 240	12	3.41	2.81 - 4.01	0.20
186395	213	174 - 252	13	3.56	2.93 - 4.19	0.21
154345	216	177 - 255	13	3.61	2.95 - 4.27	0.22
144527	246	201 - 291	15	4.11	3.36 - 4.86	0.25
132865	215	176 - 254	13	3.59	2.93 - 4.25	0.22
160386	227	185 - 269	14	3.79	3.10 - 4.48	0.23
160393	232	190 - 274	14	3.87	3.18 - 4.56	0.23
191111	239	197 - 281	14	3.99	3.27 - 4.71	0.24
191112	222	183 - 261	13	3.71	3.05 - 4.37	0.22
191125	228	186 - 270	14	3.81	3.12 - 4.50	0.23

2017-06



ALP IFCC Gen.2 (Schumann)

COBAS INTEGRA 400 plus/800

03333752 190, Sys.-ID 07 6761 1
 03333701 190, Sys.-ID 07 6760 3

ALP2S: Test-ID 0-551

ALP2S: Test-ID 0-554 (primary tube)

ALP2L: Test-ID 0-550

ALP2L: Test-ID 0-553 (primary tube)

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός
 Plasma / Plazma / Osocze / Πλάσμα

LOT	U/L	µkat/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsværdi Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

Cfas, No. 10759350 360

144452	231	3.86
144453	239	3.99
144454	239	3.99
198830	229	3.82
198831	236	3.94
198832	236	3.94
198833	236	3.94

LOT	U/L			µkat/L		
	Value/Wert	Range/Bereich	1SD/1SD	Value/Wert	Range/Bereich	1SD/1SD
	Valeur/Valor	Intervalle/Intervalo	1s/1DE	Valeur/Valor	Intervalle/Intervalo	1s/1DE
	Valore/Valor	Intervallo/Intervalo	1DS/1DP	Valore/Valor	Intervallo/Intervalo	1DS/1DP
	Värde/Værdi	Område/Område	1DP/1SD	Värde/Værdi	Område/Område	1DP/1SD
	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD	Hornota/Hodnota	Rozmezí/Rozsah	1SD/1SD
	Wartość/Érték	Zakres/Tartomány	1SD/1SD	Wartość/Érték	Zakres/Tartomány	1SD/1SD
	Τιμή	Εύρος	1SD	Τιμή	Εύρος	1SD

Precinorm U plus, No. 12149435 160

186350	85.5	70.2	-	100.8	5.1	1.43	1.16	-	1.70	0.09
203067	80.1	65.7	-	94.5	4.8	1.34	1.10	-	1.58	0.08

Precipath U plus, No. 12149443 160

186345	229	187	-	271	14	3.82	3.13	-	4.51	0.23
203060	221	182	-	260	13	3.69	3.03	-	4.35	0.22

PreciControl ClinChem Multi 1, No. 05947626 160

144514	91.7	75.2	-	108.2	5.5	1.53	1.26	-	1.80	0.09
--------	------	------	---	-------	-----	------	------	---	------	------

PreciControl ClinChem Multi 2, No. 05947774 160

144527	246	201	-	291	15	4.11	3.36	-	4.86	0.25
186395	213	174	-	252	13	3.56	2.93	-	4.19	0.21