

LOT	e-BC Ver.					U/L			µkat/L		
	c311	c501	c502	c701	c702	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

**Precinorm U, No. 10171735 122, 10171743 122**

176367	4	4	4	4	4	84.2	68.9 - 99.5	5.1	1.41	1.17 - 1.65	0.08
186436	4	4	4	4	4	89.2	73.0 - 105.4	5.4	1.49	1.22 - 1.76	0.09
186438	2	2	2	2	2	94.0	77.2 - 110.8	5.6	1.57	1.30 - 1.84	0.09
186439	2	2	2	2	2	93.4	76.6 - 110.2	5.6	1.56	1.29 - 1.83	0.09

**Precinorm U plus, No. 12149435 122**

186350	2	2	2	2	2	84.5	69.2 - 99.8	5.1	1.41	1.17 - 1.65	0.08
203067	1	1	1	1	1	75.5	62.0 - 89.0	4.5	1.26	1.02 - 1.50	0.08
203068	1	1	1	1	1	75.1	61.6 - 88.6	4.5	1.25	1.01 - 1.49	0.08

**Precipath U, No. 10171760 122, 10171778 122**

176371	4	4	4	4	4	221	182 - 260	13	3.69	3.03 - 4.35	0.22
186443	4	4	4	4	4	230	188 - 272	14	3.84	3.15 - 4.53	0.23

**Precipath U plus, No. 12149443 122**

186345	2	2	2	2	2	226	184 - 268	14	3.77	3.08 - 4.46	0.23
203060	2	2	2	2	2	208	172 - 244	12	3.47	2.84 - 4.10	0.21

**PreciControl ClinChem Multi 1, No. 05947626 190, 05117003 190, 05117208 922 (Q.C.S.)**

186373	5	5	5	5	5	98.4	80.7 - 116.1	5.9	1.64	1.34 - 1.94	0.10
186380	4	4	4	4	4	92.5	75.7 - 109.3	5.6	1.54	1.27 - 1.81	0.09
132859	2	2	2	2	2	94.5	77.4 - 111.6	5.7	1.58	1.31 - 1.85	0.09
154342	3	3	3	3	3	90.6	74.4 - 106.8	5.4	1.51	1.24 - 1.78	0.09
144514	3	3	3	3	3	89.3	73.1 - 105.5	5.4	1.49	1.22 - 1.76	0.09
160400	2	2	2	2	2	94.9	77.8 - 112.0	5.7	1.58	1.31 - 1.85	0.09
160401	2	2	2	2	2	97.6	79.9 - 115.3	5.9	1.63	1.33 - 1.93	0.10
160407	3	3	3	3	3	102	84 - 120	6	1.70	1.40 - 2.00	0.10
160410	3	3	3	3	3	102	84 - 120	6	1.70	1.40 - 2.00	0.10
191102	1	1	1	1	1	92.1	75.6 - 108.6	5.5	1.54	1.27 - 1.81	0.09
191118	1	1	1	1	1	104	86 - 122	6	1.74	1.44 - 2.04	0.10

**PreciControl ClinChem Multi 2, No. 05947774 190, 05117216 190, 05117291 922 (Q.C.S.)**

186387	5	5	5	5	5	229	187 - 271	14	3.82	3.13 - 4.51	0.23
186394	4	4	4	4	4	201	165 - 237	12	3.36	2.76 - 3.96	0.20
186395	4	4	4	4	4	213	174 - 252	13	3.56	2.93 - 4.19	0.21
154345	3	2	2	2	2	213	174 - 252	13	3.56	2.93 - 4.19	0.21
144527	3	3	3	3	3	244	199 - 289	15	4.07	3.35 - 4.79	0.24
132865	3	3	3	3	3	211	172 - 250	13	3.52	2.89 - 4.15	0.21
160386	2	2	2	2	2	222	183 - 261	13	3.71	3.05 - 4.37	0.22
160393	3	3	3	3	3	226	184 - 268	14	3.77	3.08 - 4.46	0.23
191111	1	1	1	1	1	234	192 - 276	14	3.91	3.22 - 4.60	0.23
191112	1	1	1	1	1	224	185 - 263	13	3.74	3.08 - 4.40	0.22
191125	1	1	1	1	1	225	183 - 267	14	3.76	3.07 - 4.45	0.23



### Alkaline Phosphatase acc. to IFCC Gen.2 (Schumann)

**cobas c 311/501/502 analyzers**  
 03333752 190, Sys.-ID 07 6761 1  
 03333701 190, Sys.-ID 07 6760 3

**cobas c 701/702 analyzers**  
 05166888 190, Sys.-ID 05 6760 3

For **cobas c 311/501** analyzers:  
**ALP2S:** ACN 158  
**ALP2L:** ACN 683

For **cobas c 701/702** analyzers:  
**ALP2:** ACN 8683

For **cobas c 502** analyzer:  
**ALP2S:** ACN 8158  
**ALP2L:** ACN 8683

Serum / Sérum / Suero / Siero / Soro / Sérum / Surowica / Szérum / Ορός  
 Plasma / Plazma / Osocze / Πλάσμα

LOT	U/L	µkat/L
	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärði Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης	Calibration Value/Kalibrationswert Valeur de calibration/Valor de calibración Valore di calibrazione/Valor de calibração Kalibreringsvärde/Kalibreringsvärði Kalibrační hodnota/Hodnota kalibrácie Wartość kalibratora/Kalibrációs érték Τιμή βαθμονόμησης

Cfas, No. 10759350 360

144452	233	3.89
144453	239	3.99
144454	239	3.99
198830	225	3.76
198831	233	3.89
198832	233	3.89
198833	233	3.89

# ALP2

cobas c 311/501/502/701/702 analyzers

Note 2017-06 - 03333752 - 2 - 4

LOT	U/L			µkat/L		
	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD	Value/Wert Valeur/Valor Valore/Valor Värde/Værdi Hornota/Hodnota Wartość/Érték Τιμή	Range/Bereich Intervalle/Intervalo Intervallo/Intervalo Område/Område Rozmezí/Rozsah Zakres/Tartomány Εύρος	1SD/1SD 1s/1DE 1DS/1DP 1DP/1SD 1SD/1SD 1SD/1SD 1SD

### Precinorm U plus, No. 12149435 160

186350	84.5	69.2 - 99.8	5.1	1.41	1.17 - 1.65	0.08
203067	75.5	62.0 - 89.0	4.5	1.26	1.02 - 1.50	0.08

### Precipath U plus, No. 12149443 160

186345	226	184 - 268	14	3.77	3.08 - 4.46	0.23
203060	208	172 - 244	12	3.47	2.84 - 4.10	0.21

### PreciControl ClinChem Multi 1, No. 05947626 160

144514	89.3	73.1 - 105.5	5.4	1.49	1.22 - 1.76	0.09
--------	------	--------------	-----	------	-------------	------

### PreciControl ClinChem Multi 2, No. 05947774 160

144527	244	199 - 289	15	4.07	3.35 - 4.79	0.24
186395	213	174 - 252	13	3.56	2.93 - 4.19	0.21